

AAAS

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 59kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human AAAS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AAA; AAASb; GL003; ALADIN; ADRACALA; ADRACALIN;

Catalog #: A6427

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8086

Isotype: IgG

Swiss Prot: Q9NRG9

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the WD-repeat family of regulatory proteins and may be involved in normal development of the peripheral and central nervous system. The encoded protein is part of the nuclear pore complex and is anchored there by NDC1. Defects in this gene are a cause of achalasia-addisonianism-alacrima syndrome (AAAS), also called triple-A syndrome or Allgrove syndrome. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).