

C21orf33

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 28kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human C21orf33

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ES1; HES1; KNPH; KNPI; GT335;

Catalog #: A6429

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8209

Isotype: IgG

Swiss Prot: P30042

Purity: Affinity purification

For research use only.

Background:

This gene encodes a potential mitochondrial protein that is a member of the DJ-1/Pfpl gene family.

This protein is overexpressed in fetal Down syndrome brain. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).