

NAPG

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 34kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NAPG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GAMMASNAP;

Catalog #: A6432

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8774

Isotype: IgG

Swiss Prot: Q99747

Purity: Affinity purification

For research use only.

Background:

This gene encodes soluble NSF attachment protein gamma. The soluble NSF attachment proteins (SNAPs) enable N-ethyl-maleimide-sensitive fusion protein (NSF) to bind to target membranes.

NSF and SNAPs appear to be general components of the intracellular membrane fusion apparatus, and their action at specific sites of fusion must be controlled by SNAP receptors particular to the membranes being fused. The product of this gene mediates platelet exocytosis and controls the membrane fusion events of this process.

To place an order, please [Click HERE](#).