

CD84

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 38kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CD84

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LY9B; hCD84; mCD84; SLAMF5;

Catalog #: A6433

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8832

Isotype: IgG

Swiss Prot: Q9UIB8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a membrane glycoprotein that is a member of the signaling lymphocyte activation molecule (SLAM) family. This family forms a subset of the larger CD2 cell-surface receptor Ig superfamily. The encoded protein is a homophilic adhesion molecule that is expressed in numerous immune cells types and is involved in regulating receptor-mediated signaling in those cells. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).