

SKAP2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 41kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SKAP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PRAP; RA70; SAPS; SCAP2; SKAP55R; SKAP-HOM;

Catalog #: A6435

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8935

Isotype: IgG

Swiss Prot: O75563

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the src family kinases. This protein is similar to the src kinase associated phosphoprotein 1. It is an adaptor protein that is thought to play an essential role in the src signaling pathway in various cells. It inhibits PTK2B/RAFTK activity and regulates alpha-synuclein phosphorylation.

To place an order, please [Click HERE](#).