

## FIBP

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 41kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human FIBP

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

FGFIBP; FIBP-1;

**Catalog #:** A6436

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 9158

**Isotype:** IgG

**Swiss Prot:** O43427

**Purity:** Affinity purification

For research use only.

**Background:**

Acidic fibroblast growth factor is mitogenic for a variety of different cell types and acts by stimulating mitogenesis or inducing morphological changes and differentiation. The FIBP protein is an intracellular protein that binds selectively to acidic fibroblast growth factor (aFGF). It is postulated that FIBP may be involved in the mitogenic action of aFGF. Two transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*