

DHRS2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 29kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human DHRS2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HEP27; SDR25C1;

Catalog #: A6446

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10202

Isotype: IgG

Swiss Prot: Q13268

Purity: Affinity purification

For research use only.

Background:

DHRS2 (dehydrogenase/reductase (SDR family) member 2) is a protein-coding gene. Diseases associated with DHRS2 include otitis media, and hepatoblastoma. GO annotations related to this gene include carbonyl reductase (NADPH) activity. An important paralog of this gene is HPGD.

To place an order, please [Click HERE](#).