

LYVE1

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 35kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LYVE1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HAR; XLKD1; LYVE-1; CRSBP-1;

Catalog #: A6452

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10894

Isotype: IgG

Swiss Prot: Q9Y5Y7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a type I integral membrane glycoprotein. The encoded protein acts as a receptor and binds to both soluble and immobilized hyaluronan. This protein may function in lymphatic hyaluronan transport and have a role in tumor metastasis.

To place an order, please [Click HERE](#).