

SEPHS1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 42kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SEPHS1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SPS; SELD; SPS1;

Catalog #: A6454

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 22929

Isotype: IgG

Swiss Prot: P49903

Purity: Affinity purification

For research use only.

Background:

This gene encodes an enzyme that synthesizes selenophosphate from selenide and ATP. Selenophosphate is the selenium donor used to synthesize selenocysteine, which is co-translationally incorporated into selenoproteins at in-frame UGA codons.

To place an order, please [Click HERE](#).