

DDAH2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 30kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human DDAH2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

G6a; DDAH; NG30; DDAHII; HEL-S-277;

Catalog #: A6457

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23564

Isotype: IgG

Swiss Prot: O95865

Purity: Affinity purification

For research use only.

Background:

This gene belongs to the dimethylarginine dimethylaminohydrolase (DDAH) gene family. The encoded enzyme plays a role in nitric oxide generation by regulating cellular concentrations of methylarginines, which in turn inhibit nitric oxide synthase activity.

To place an order, please [Click HERE](#).