

GNL3

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 61kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GNL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NS; E2IG3; NNP47; C77032;

Catalog #: A6459

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26354

Isotype: IgG

Swiss Prot: Q9BVP2

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene may interact with p53 and may be involved in tumorigenesis.

The encoded protein also appears to be important for stem cell proliferation. This protein is found in both the nucleus and nucleolus. Three transcript variants encoding two different isoforms have been found for this gene.

To place an order, please [Click HERE](#).