

NRBF2

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 32kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NRBF2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

COPR1; COPR2; NRBF-2;

To place an order, please [Click HERE](#).

Catalog #: A6462

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29982

Isotype: IgG

Swiss Prot: Q96F24

Purity: Affinity purification

For research use only.