

POLE3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 16kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human POLE3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

p17; YBL1; CHRAC17; CHARAC17;

Catalog #: A6469

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54107

Isotype: IgG

Swiss Prot: Q9NRF9

Purity: Affinity purification

For research use only.

Background:

POLE3 is a histone-fold protein that interacts with other histone-fold proteins to bind DNA in a sequence-independent manner. These histone-fold protein dimers combine within larger enzymatic complexes for DNA transcription, replication, and packaging.

To place an order, please [Click HERE](#).