

KRBOX4

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 20kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human KRBOX4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ZNF673;

Catalog #: A6471

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55634

Isotype: IgG

Swiss Prot: Q5JUW0

Purity: Affinity purification

For research use only.

Background:

This encodes a zinc finger protein with an N-terminal KRAB (Kruppel-associated) domain found in transcriptional repressors. This gene is located in a region of the X chromosome thought to be involved in nonsyndromic X-linked mental retardation. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).