

SMYD2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 49kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SMYD2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

KMT3C; HSKM-B; ZMYND14;

Catalog #: A6474

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 56950

Isotype: IgG

Swiss Prot: Q9NRG4

Purity: Affinity purification

For research use only.

Background:

SET domain-containing proteins, such as SMYD2, catalyze lysine methylation (Brown et al., 2006

To place an order, please [Click HERE](#).