

CYP20A1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 52kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CYP20A1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CYP-M;

Catalog #: A6476

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57404

Isotype: IgG

Swiss Prot: Q6UW02

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases that catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein lacks one amino acid of the conserved heme binding site. It also lacks the conserved I-helix motif AGX(D,E)T, suggesting that its substrate may carry its own oxygen.

To place an order, please [Click HERE](#).