

NGB

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 16kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NGB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6477

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 58157

Isotype: IgG

Swiss Prot: Q9NPG2

Purity: Affinity purification

For research use only.

Background:

This gene encodes an oxygen-binding protein that is distantly related to members of the globin gene family. It is highly conserved among other vertebrates. It is expressed in the central and peripheral nervous system where it may be involved in increasing oxygen availability and providing protection under hypoxic/ischemic conditions.

To place an order, please [Click HERE](#).