

BCL2L14

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 36kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human BCL2L14

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BCLG;

Catalog #: A6478

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79370

Isotype: IgG

Swiss Prot: Q9BZR8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the BCL2 protein family. BCL2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. Overexpression of this gene has been shown to induce apoptosis in cells. Three alternatively spliced transcript variants encoding two distinct isoforms have been reported for this gene.

To place an order, please [Click HERE](#).