

KIF2B

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 76kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human KIF2B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6480

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84643

Isotype: IgG

Swiss Prot: Q8N4N8

Purity: Affinity purification

For research use only.

Background:

KIF2B (kinesin family member 2B) is a protein-coding gene. Diseases associated with KIF2B include charcot-marie-tooth disease type 4, and hereditary neuropathies. GO annotations related to this gene include microtubule binding and microtubule motor activity. An important paralog of this gene is KIF2A.

To place an order, please [Click HERE](#).