

LNX1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 80kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human LNX1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LNX; MPDZ; PDZRN2;

Catalog #: A6481

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84708

Isotype: IgG

Swiss Prot: Q8TBB1

Purity: Affinity purification

For research use only.

Background:

This gene encodes a membrane-bound protein that is involved in signal transduction and protein interactions. The encoded product is an E3 ubiquitin-protein ligase, which mediates ubiquitination and subsequent proteasomal degradation of proteins containing phosphotyrosine binding (PTB) domains. This protein may play an important role in tumorigenesis. Alternatively spliced transcript variants encoding distinct isoforms have been described. A pseudogene, which is located on chromosome 17, has been identified for this gene.

To place an order, please [Click HERE](#).