

POMGNT2

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50

Calculated MW: 66kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human POMGNT2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AGO61; GTDC2; MDDGA8; C3orf39;

Catalog #: A6484

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84892

Isotype: IgG

Swiss Prot: Q8NAT1

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein with glycosyltransferase activity although its function is not currently known.

To place an order, please [Click HERE](#).