

RFFL

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 40kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human RFFL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CARP2; FRING; CARP-2; RNF189; RNF34L; RIFIFYLIN;

To place an order, please [Click HERE](#).

Catalog #: A6489

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 117584

Isotype: IgG

Swiss Prot: Q8WZ73

Purity: Affinity purification

For research use only.