

FLCN

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 64kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human FLCN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BHD; FLCL;

Catalog #: A6493

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 201163

Isotype: IgG

Swiss Prot: Q8NFG4

Purity: Affinity purification

For research use only.

Background:

This gene is located within the Smith-Magenis syndrome region on chromosome 17. Mutations in this gene are associated with Birt-Hogg-Dube syndrome, which is characterized by fibrofolliculomas, renal tumors, lung cysts, and pneumothorax. Alternative splicing of this gene results in two transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).