

ZBTB38

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 134kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ZBTB38

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CIBZ; ZNF921; PPP1R171;

Catalog #: A6494

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 253461

Isotype: IgG

Swiss Prot: Q8NAP3

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a zinc finger transcriptional activator that binds methylated DNA. The encoded protein can form homodimers or heterodimers through the zinc finger domains. In mouse, inhibition of this protein has been associated with apoptosis in some cell types.

To place an order, please [Click HERE](#).