

ARL3

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 20kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic Peptide of human ARL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

ARL3;

Catalog #: A6499

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 403

Isotype: IgG

Swiss Prot: P36405

Purity: Affinity purification

For research use only.

Background:

ADP-ribosylation factor-like 3 is a member of the ADP-ribosylation factor family of GTP-binding proteins. ARL3 binds guanine nucleotides but lacks ADP-ribosylation factor activity.

To place an order, please [Click HERE](#).