

ASF1B

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 22kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ASF1B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CIA-II;

Catalog #: A6511

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55723

Isotype: IgG

Swiss Prot: Q9NVP2

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the H3/H4 family of histone chaperone proteins and is similar to the anti-silencing function-1 gene in yeast. The encoded protein is the substrate of the tousel-like kinase family of cell cycle-regulated kinases, and may play a key role in modulating the nucleosome structure of chromatin by ensuring a constant supply of histones at sites of nucleosome assembly.

To place an order, please [Click HERE](#).