

ADSS

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 50kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ADSS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ADEH; ADSS 2;

Catalog #: A6516

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 159

Isotype: IgG

Swiss Prot: P30520

Purity: Affinity purification

For research use only.

Background:

This gene encodes the enzyme adenylosuccinate synthetase which catalyzes the first committed step in the conversion of inosine monophosphate to adenosine monophosphate. A pseudogene of this gene is found on chromosome 17.

To place an order, please [Click HERE](#).