

ANLN

Reactivity: Human

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW: 124kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ANLN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

scra; FSGS8; Scraps;

Catalog #: A6524

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54443

Isotype: IgG

Swiss Prot: Q9NQW6

Purity: Affinity purification

For research use only.

Background:

This gene encodes an actin-binding protein that plays a role in cell growth and migration, and in cytokinesis. The encoded protein is thought to regulate actin cytoskeletal dynamics in podocytes, components of the glomerulus. Mutations in this gene are associated with focal segmental glomerulosclerosis 8. Alternative splicing results in multiple transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).