

ANTXR2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 54kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ANTXR2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HFS; ISH; JHF; CMG2; CMG-2;

Catalog #: A6526

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 118429

Isotype: IgG

Swiss Prot: P58335

Purity: Affinity purification

For research use only.

Background:

This gene encodes a receptor for anthrax toxin. The protein binds to collagen IV and laminin, suggesting that it may be involved in extracellular matrix adhesion. Mutations in this gene cause juvenile hyaline fibromatosis and infantile systemic hyalinosis. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).