

ASF1A

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 22kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ASF1A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CIA; CGI-98; HSPC146;

Catalog #: A6528

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 25842

Isotype: IgG

Swiss Prot: Q9Y294

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the H3/H4 family of histone chaperone proteins and is similar to the anti-silencing function-1 gene in yeast. The protein is a key component of a histone donor complex that functions in nucleosome assembly. It interacts with histones H3 and H4, and functions together with a chromatin assembly factor during DNA replication and repair.

To place an order, please [Click HERE](#).