

ATP6AP2

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ATP6AP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PRR; M8-9; MRXE; RENR; XMRE; XPDS; HT028; MRXSH; ELDF10; ATP6IP2; MSTP009;
APT6M8-9; ATP6M8-9;

Catalog #: A6531

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10159

Isotype: IgG

Swiss Prot: O75787

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that is associated with adenosine triphosphatases (ATPases). Proton-translocating ATPases have fundamental roles in energy conservation, secondary active transport, acidification of intracellular compartments, and cellular pH homeostasis. There are three classes of ATPases- F, P, and V. The vacuolar (V-type) ATPases have a transmembrane proton-conducting sector and an extramembrane catalytic sector. The encoded protein has been found associated with the transmembrane sector of the V-type ATPases.

To place an order, please [Click HERE](#).