

INTS10

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 82kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human INTS10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

INT10; C8orf35;

Catalog #: A6632

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55174

Isotype: IgG

Swiss Prot: Q9NVR2

Purity: Affinity purification

For research use only.

Background:

INTS10 is a subunit of the Integrator complex, which associates with the C-terminal domain of RNA polymerase II large subunit (POLR2A; MIM 180660) and mediates 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM 180680) and U2 (RNU2; MIM 180690) (Baillat et al., 2005)

To place an order, please [Click HERE](#).