

IRAK2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:200

Calculated MW: 69kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IRAK2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IRAK-2;

Catalog #: A6635

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3656

Isotype: IgG

Swiss Prot: O43187

Purity: Affinity purification

For research use only.

Background:

IRAK2 encodes the interleukin-1 receptor-associated kinase 2, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. IRAK2 is reported to participate in the IL1-induced upregulation of NF-kappaB.

To place an order, please [Click HERE](#).