

KHDRBS3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human KHDRBS3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Etle; SALP; SLM2; SLM-2; TSTAR; T-STAR; etoile;

To place an order, please [Click HERE](#).

Catalog #: A6637

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10656

Isotype: IgG

Swiss Prot: O75525

Purity: Affinity purification

For research use only.