

KIF1B

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 204kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human KIF1B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KLP; CMT2; CMT2A; CMT2A1; HMSNII; NBLST1;

Background:

This gene encodes a motor protein that transports mitochondria and synaptic vesicle precursors. Mutations in this gene cause Charcot-Marie-Tooth disease, type 2A1.

To place an order, please [Click HERE](#).

Catalog #: A6638

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23095

Isotype: IgG

Swiss Prot: O60333

Purity: Affinity purification

For research use only.