

LGALS8

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 35kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LGALS8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Gal-8; PCTA1; PCTA-1; Po66-CBP;

Catalog #: A6646

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3964

Isotype: IgG

Swiss Prot: O00214

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the galectin family. Galectins are beta-galactoside-binding animal lectins with conserved carbohydrate recognition domains. The galectins have been implicated in many essential functions including development, differentiation, cell-cell adhesion, cell-matrix interaction, growth regulation, apoptosis, and RNA splicing. This gene is widely expressed in tumoral tissues and seems to be involved in integrin-like cell interactions. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).