

## MFGE8

---

**Reactivity:** Human

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 43kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human MFGE8

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

BA46; HMFG; MFGM; SED1; hP47; EDIL1; MFG-E8; SPAG10; OAcGD3S; HsT19888;

*To place an order, please [Click HERE](#).*

**Catalog #:** A6652

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 4240

**Isotype:** IgG

**Swiss Prot:** Q08431

**Purity:** Affinity purification

For research use only.