

MF12

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:80kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human MF12

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MTF1; CD228; MAP97;

Catalog #:A6653

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4241

Isotype:IgG

Swiss Prot:P08582

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a cell-surface glycoprotein found on melanoma cells. The protein shares sequence similarity and iron-binding properties with members of the transferrin superfamily. The importance of the iron binding function has not yet been identified. This gene resides in the same region of chromosome 3 as members of the transferrin superfamily. Alternative splicing results in two transcript variants.

To place an order, please [Click HERE](#).