

MINA

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 52kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MINA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

ROX; MDIG; NO52; MINA53;

Catalog #: A6655

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84864

Isotype: IgG

Swiss Prot: Q8IUJ8

Purity: Affinity purification

For research use only.

Background:

MINA is a c-Myc (MYC; MIM 190080) target gene that may play a role in cell proliferation or regulation of cell growth. (Tsuneoka et al., 2002)

To place an order, please [Click HERE](#).