

MMRN1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 138kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MMRN1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ECM; MMRN; GPIa*; EMILIN4;

Catalog #: A6658

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 22915

Isotype: IgG

Swiss Prot: Q13201

Purity: Affinity purification

For research use only.

Background:

Multimerin is a massive, soluble protein found in platelets and in the endothelium of blood vessels.

It is comprised of subunits linked by interchain disulfide bonds to form large, variably sized homomultimers. Multimerin is a factor V/Va-binding protein and may function as a carrier protein for platelet factor V. It may also have functions as an extracellular matrix or adhesive protein.

Recently, patients with an unusual autosomal-dominant bleeding disorder (factor V Quebec) were found to have a deficiency of platelet multimerin.

To place an order, please [Click HERE](#).