

MPP2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 64kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MPP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DLG2;

Catalog #: A6659

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4355

Isotype: IgG

Swiss Prot: Q14168

Purity: Affinity purification

For research use only.

Background:

Palmitoylated membrane protein 2 is a member of a family of membrane-associated proteins termed MAGUKs (membrane-associated guanylate kinase homologs). MAGUKs interact with the cytoskeleton and regulate cell proliferation, signaling pathways, and intracellular junctions.

Palmitoylated membrane protein 2 contains a conserved sequence, called the SH3 (src homology 3) motif, found in several other proteins that associate with the cytoskeleton and are suspected to play important roles in signal transduction.

To place an order, please [Click HERE](#).