

NFATC3

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 115kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NFATC3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NFAT4; NFATX;

Catalog #: A6666

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4775

Isotype: IgG

Swiss Prot: Q12968

Purity: Affinity purification

For research use only.

Background:

The product of this gene is a member of the nuclear factors of activated T cells DNA-binding transcription complex. This complex consists of at least two components: a preexisting cytosolic component that translocates to the nucleus upon T cell receptor (TCR) stimulation and an inducible nuclear component. Other members of this family participate to form this complex also. The product of this gene plays a role in the regulation of gene expression in T cells and immature thymocytes. Several transcript variants encoding distinct isoforms have been identified for this gene.

To place an order, please [Click HERE](#).