

NR4A1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 64kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NR4A1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HMR; N10; TR3; NP10; GFRP1; NAK-1; NGFIB; NUR77;

Catalog #: A6676

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3164

Isotype: IgG

Swiss Prot: P22736

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the steroid-thyroid hormone-retinoid receptor superfamily. Expression is induced by phytohemagglutinin in human lymphocytes and by serum stimulation of arrested fibroblasts. The encoded protein acts as a nuclear transcription factor. Translocation of the protein from the nucleus to mitochondria induces apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).