

NUDC

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 38kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NUDC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HNUDC; MNUDC; NPD011;

Catalog #: A6678

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10726

Isotype: IgG

Swiss Prot: Q9Y266

Purity: Affinity purification

For research use only.

Background:

This gene encodes a nuclear distribution protein that plays an essential role in mitosis and cytokinesis. The encoded protein is involved in spindle formation during mitosis and in microtubule organization during cytokinesis. Pseudogenes of this gene are found on chromosome 2.

To place an order, please [Click HERE](#).