

OGN

Reactivity: Human Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human OGN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OG; OIF; SLRR3A;

Catalog #: A6679

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4969

Isotype: IgG

Swiss Prot: P20774

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein which induces ectopic bone formation in conjunction with transforming growth factor beta. This protein is a small proteoglycan which contains tandem leucine-rich repeats (LRR). Alternatively spliced transcript variants that encode the same protein have been found for this gene.

To place an order, please [Click HERE](#).