

PRPF19

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 55kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PRPF19

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

PSO4; SNEV; PRP19; UBOX4; hPSO4; NMP200;

Catalog #: A6706

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27339

Isotype: IgG

Swiss Prot: Q9UMS4

Purity: Affinity purification

For research use only.

Background:

PSO4 is the human homolog of yeast Pso4, a gene essential for cell survival and DNA repair (Beck et al., 2008)

To place an order, please [Click HERE](#).