

QPCT

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 40kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human QPCT

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

QC; GCT; sQC;

Catalog #: A6711

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 25797

Isotype: IgG

Swiss Prot: Q16769

Purity: Affinity purification

For research use only.

Background:

This gene encodes human pituitary glutaminy cyclase, which is responsible for the presence of pyroglutamyl residues in many neuroendocrine peptides. The amino acid sequence of this enzyme is 86% identical to that of bovine glutaminy cyclase.

To place an order, please [Click HERE](#).