

RGS3

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 132kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RGS3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C2PA; RGP3;

Catalog #: A6719

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5998

Isotype: IgG

Swiss Prot: P49796

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the regulator of G-protein signaling (RGS) family. This protein is a GTPase-activating protein that inhibits G-protein-mediated signal transduction. Alternative splicing and the use of alternative promoters results in multiple transcript variants encoding different isoforms. Long isoforms are largely cytosolic and plasma membrane-associated with a function in Wnt signaling and in the epithelial mesenchymal transition, while shorter N-terminally-truncated isoforms can be nuclear.

To place an order, please [Click HERE](#).