

RPH3A

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 76kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RPH3A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6722

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 22895

Isotype: IgG

Swiss Prot: Q9Y2J0

Purity: Affinity purification

For research use only.

Background:

Exocytosis of neurotransmitters and hormones is fundamental to synaptic neurotransmission and cell-cell communication. RAB3A (MIM 179390) is a small G protein that is thought to act at late stages of exocytosis, and RPH3A is a RAB3A effector (Lin et al., 2007)

To place an order, please [Click HERE](#).