

SPINT2

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 28kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SPINT2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PB; Kop; HAI2; DIAR3; HAI-2;

Catalog #: A6749

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10653

Isotype: IgG

Swiss Prot: O43291

Purity: Affinity purification

For research use only.

Background:

This gene encodes a transmembrane protein with two extracellular Kunitz domains that inhibits a variety of serine proteases. The protein inhibits HGF activator which prevents the formation of active hepatocyte growth factor. This gene is a putative tumor suppressor, and mutations in this gene result in congenital sodium diarrhea. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).